



# Strategies to Evaluate CCI of Vial and Syringe Systems over Time and Temperature

**Amy Gindhart, Louis Brasten, Amy Kim, Jennifer Roark**

West Pharmaceutical Services, Inc.






Proper characterization of container closure integrity (CCI) of vial and syringe systems is essential



Aid pharmaceutical manufacturers to determine if a system is suitable for a drug product



Determine if a system can meet the maximum allowable leakage limit (MALL)



No single method is suitable for all systems



Proper characterization strategy should employ multiple methods



Long-term case studies using tracer gas detection vacuum mode (helium leak detection) and laser-based gas headspace analysis (oxygen and carbon dioxide)



Glass/elastomer vial combinations over 2 years at room temperature



Polymer/elastomer vial combinations over 2 years at room temperature, ultra-cold, and cryogenic, including carbon dioxide atmosphere



Glass syringe combinations



Discussed also are the specific utilities of each method



Techniques combined offer more complete picture of CCI and aid in determining if MALL is met



## Residual Seal Force (complementary method)

- Quantifies compression of stopper against vial



## Tracer Gas detection vacuum mode with Helium (He Leak)

- Measures stopper/vial seal integrity
- Deterministic – endorsed by USP <1207>



## Laser-Based Gas Headspace Analysis (O<sub>2</sub> & CO<sub>2</sub>)

- Measures diffusion/effusion through vial, stopper, and interface
- Deterministic – endorsed by USP <1207>

# Residual Seal Force (RSF)



Equipment: Genesis Model AWG RSF Tester



RSF is a measurement of the seal tightness of the stopper against the vial in the stopper/seal/vial combination resulting from the sealing process



Indirectly measures the force exerted by the stopper on a vial's land surface



Sealed vial is placed into the instrument holder with an appropriately sized cap anvil on top of the vial




Pressure is applied to the sample; the machine calculates the force (pounds-force) required to dislodge the crimp seal from the underside of the vial crown.



# Helium Leak

mass spectrometer



Room Temperature, Ultra-low & Cryogenic Temperatures



Equipment: Leak Detection Associates (LDA)/Packaging Technologies and Inspection (PTI) 1284+ Seal Integrity Monitoring System (SIMS)



2 modes of operation

- Standard (room temperature,  $-80^{\circ}\text{C}$  and  $-180^{\circ}\text{C}$ )
- Continuous flow (room temperature)




Mass spec measures level of Helium and reports as leak rate (mbar-L/s)



Advantages

- Highly sensitive
- Quantitative leak size determination
- Non-destructive
- 100% testing feasible

mass spectrometer




Room Temperature



Often used as a valuation standard for CCI performance

Hole Diameter ( $\mu$ )	He-Leak ( $\text{cm}^3/\text{s}$ @ STP)	- log (He-Leak)	Microbial Ingress Rate (%)
2	$1.0 \times 10^{-3}$	3.0	70
0.7	$2.0 \times 10^{-4}$	3.7	65
0.4	$9.0 \times 10^{-6}$	5.0	11
---	$6 \times 10^{-6}$	5.2	8 - 10
0.3	$2.0 \times 10^{-6}$	5.7	7
0.2	$2.2 \times 10^{-7}$	6.6	0
0.1	$1.0 \times 10^{-7}$	7.0	0



<sup>1</sup> L. E. Kirsch, et al. (University of Iowa) *Pharmaceutical Container/Closure Integrity II: The Relationship between Microbial Ingress and Helium Leak Rates in Rubber-Stoppered Vials*. PDA Journal of Pharmaceutical Science & Technology, 51 (5), 195-202 (1997).

# Headspace Analysis



Equipment: Oxygen  
Headspace Analyzer, FMS-760  
by Lighthouse Instruments



Amount absorbed correlates  
to concentration





Uses near-IR radiation  
wavelength that is absorbed  
by gas of interest



Advantages

- Rapid testing time
- Non-Destructive
- Highly sensitive
- No sample preparation



# Glass Vial / Stopper / Seal Combinations



He leak and O<sub>2</sub> headspace routinely used at West to quantify performance of vial/stopper/seal combinations varying in:



Vial: size, blowback, supplier



Stopper: size, configuration/design, elastomer, post-treatment



Assembly: compression level, seal



Example Presented:




13 mm serum bromobutyl elastomer stopper



2R straightwall glass vial



# Residual Seal Force at Ambient Temperature



2R straight wall glass vial with 13mm serum bromobutyl elastomer stopper




Increases with compression level as expected



Reduces in value quickly at 3M and then remains constant over 24M due to the relaxation of the elastomer

# Helium Leak at Ambient Temperature



There was no change in helium leak values observed over 24M



Data demonstrates excellent sealing integrity




Values correspond to 0% microbial ingress per Kirsch limit

# Helium Leak at Ambient Temperature (Low Compression)





Helium leak can indicate compression levels that are not optimal

- Including standard deviation



# Headspace Analysis (O<sub>2</sub>) at Ambient Temperature



There was some egress of O<sub>2</sub> and ingress of N<sub>2</sub> over time not predicted by the helium leak data.



Rubber is known to be gas permeable



Egress was on avg 1.7% O<sub>2</sub> after 24M.

# Conclusion of Glass/Elastomer Vial Combinations




Combined He leak and O<sub>2</sub> headspace analyses give a more complete view of CCI and better enables understanding if MALL can be met




For the present example (13 mm stopper and 2R vial)

- Seal integrity over 2 years is excellent
- All He leak values are better than the Kirsch limit (correspond to ~ 0% microbial ingress)
- O<sub>2</sub> headspace values show some quantified permeation through the rubber over time



# Reason for Low Temp Interest in Polymer Vials



Gene and cell therapies require storage at  $-80^{\circ}\text{C}$  (ultra-low) and below  $-130^{\circ}\text{C}$  (cryogenic).




The vial/stopper/seal combinations selected must be able to maintain container closure integrity (CCI) during long-term storage at said temperatures.



Polymer-based vial/stopper/seal combinations have lower risk of loss of CCI at low temperature


- Based on better match of Coefficient of Thermal Expansion (CTE)

# Polymer Vial-Stopper-Seal Combinations





**2005 Nobel Prize**

Ring-Opening Metathesis Polymerization  
e.g., Ru[P<sub>2</sub>C<sub>3</sub>H<sub>2</sub>(C<sub>6</sub>H<sub>5</sub>)<sub>8</sub>]Cl<sub>2</sub>




Hydrogenation



Mw ~ 7 x 10<sup>4</sup>  
PD ~ 2

Daikyo Crystal Zenith® Cyclic Olefin Polymer (COP)

# Differences in Coefficient of Thermal Expansion






# Differences in Coefficient of Thermal Expansion



Room  
Temperature



Stopper



Vial

	Coefficient of Thermal Expansion <sup>1</sup> ( $10^{-6}$ cm / cm – K)	Volume Shrinkage (%) (room temperature to -180°C)
Glass	4	0.2
Rubber	77	4.5
COP	70 <sup>2</sup>	4.1

<sup>1</sup>Thermal Expansion-Linear Expansion Coefficients: [http://www.engineeringtoolbox.com/linear-expansion-coefficients-d\\_95.html](http://www.engineeringtoolbox.com/linear-expansion-coefficients-d_95.html)

<sup>2</sup>Characteristics of Daikyo Resin CZ, Daikyo Seiko, Ltd. Technical Report DS-CZ-E017, January 2022.

# Oxygen Headspace Room Temperature



2 ml COP vial and 13 mm serum stopper



Significant permeation of O<sub>2</sub>




Expected –polymers are gas permeable

# Oxygen Headspace -80°C Storage



Oxygen Headspace for -80°C Storage (N<sub>2</sub> filled vials)



2 ml COP vial and 13 mm serum stopper




Essentially no ingress of O<sub>2</sub> over 12 M



Low level of ingress of O<sub>2</sub> over 24 M

# Oxygen Headspace $\leq -130^{\circ}\text{C}$ Storage



2 mL COP vial and 13 mm serum stopper




No ingress of  $\text{N}_2$  over 24 M

# CO<sub>2</sub> Headspace at -78°C (Dry Ice Storage)



2R COP vial and 13 mm serum stopper



Upon storage on dry ice and thaw small amount of CO<sub>2</sub> appears – results from dissolution/permeation of CO<sub>2</sub> in COP



Level is low and quantified



Prevented by secondary packaging (poly(ethylene terephthalate)-based sealed bag)

# Conclusion of Polymer/Elastomer Vial Combinations



For COP-based systems, gas permeation is quantified:

- Significant at room temperature
- Substantially reduced at  $-80^{\circ}\text{C}$
- Eliminated at  $\leq -130^{\circ}\text{C}$



For COP-based systems stored on dry ice:


- Small level of  $\text{CO}_2$  ingress occurs
- Can be addressed with a secondary package





Combining techniques (RSF, tracer gas analysis and headspace analysis) enables a more complete and quantified measure of CCI for vial and syringe combinations to enable judgement if MALL can be met

---



## Project Managers & Analytical Laboratory Staff:

- V. Gupta, B. Jacobs, E. Crouch, G. O'deens, A. Chapman, A. McLean

## References:

- L. E. Kirsch, et al. (University of Iowa) *Pharmaceutical Container/Closure Integrity II: The Relationship between Microbial Ingress and Helium Leak Rates in Rubber-Stoppered Vials*. PDA Journal of Pharmaceutical Science & Technology, 51 (5), 195-202 (1997).
- Gehrman M.R., McAndrew T.P. Long-Term Study of Container Closure Integrity of Rubber-Glass Vial Systems by Multiple Methods. PDA J Pharm Sci Technol. 2020 Jan-Feb;74(1):147-161. doi: 10.5731/pdajpst.2019.010223
- J. Edwards, W. Wolf, R. Grubbs. *The synthesis of cyclic polymers by olefin metathesis: Achievements and challenges*. J. Polym. Sci., Part A: Polym. Chem. 2019, 57, 228–242.
- Ingress of Gases into Cyclic Olefin Polymer Vial-Based Container Closure Systems at -80°C, West Pharmaceuticals Service Inc. Technical Report 2019/028
- Long-Term Container Closure Integrity Testing of Vial-Stopper-Seal Combinations Comprising Daikyo Crystal Zenith® Vials at Ultra-Low and Cryogenic Temperatures, West Pharmaceuticals Service Inc. Technical Report 2023/261
- M. Gehrman, O. Laskina, L. Fang, P. McAndrew. Challenges in Low-Temperature Storage of Cell Therapy Drug Products. PepTalk 2020, San Diego, CA, January 20, 2020.
- Characteristics of Daikyo Resin CZ, Daikyo Seiko, Ltd. Technical Report DS-CZ-E017, January 2022.