



Elevate **standardization**
and **sterility**.

adaptiQ® ready-to-use pharmaceutical glass vials

Discover standardized flexibility with SCHOTT Pharma's adaptiQ® – delivering gold standard borosilicate glass vials in a pre-washed, pre-sterilized, ready-to-use (RTU) format. Part of the SCHOTT iQ® platform, the global RTU standard, adaptiQ® is packaged in the industry standard 3" tub and nest configuration. By eliminating glass-to-glass contact, adaptiQ® prevents breakage, maintains high quality, and lowers the total cost of ownership (TCO).


Developed in collaboration with leading machine suppliers, adaptiQ® seamlessly integrates into various new and existing filling lines. Its standardization, based on the decade-proven ISO 11040-7 and other relevant industry standards such as ISO 21882, significantly reduces changeover times.

Maximize flexibility with adaptiQ®

- Easy access to high-quality vials for clinical and commercial fills
- Seamless switch between container formats, allowing different vial types and formats to be filled on a single line
- Integrated advanced washing and sterilization processes for high-quality vials

Increase standardization with adaptiQ®

- Compatible with all SCHOTT Pharma bulk vials for maximum flexibility
- Compliant with all relevant industry standards for reduced complexity
- Minimized particle load, glass-to-glass contact, and container breakage increases patient safety
- High-quality packaging materials result in an exceptionally reliable solution



Upgrade
your filling
operation today.

Discover more at:
schott-pharma.com

schott-pharma.com

SCHOTT Pharma AG & Co. KGaA, Hattenbergstraße 10, 55122 Mainz, Germany
pharma@schott.com

SCHOTT
PHARMA