



MixJect® Transfer Device

Integrated device for drug transfer and injection



The MixJect® Transfer Device enables easy preparation of lyophilized drugs. The transfer device is intended for the transfer and injection of drugs contained in a vial. Upon transfer, the drug is available for immediate injection with a pre-attached needle.

MIXJECT® TRANSFER DEVICE ATTRIBUTE

- Preassembled with a client-specified standard safety needle
- Available in two sizes: 13mm and 20mm
- Fits standard vials, syringes, needles having Luer lock connections

MIXJECT® TRANSFER DEVICE IS IDEAL FOR

- Lyophilized drugs supplied with prefilled diluent syringes

INNOVATIONS IN SAFETY AND ADMINISTRATION SYSTEMS

West's administration systems enhance the safe delivery of drugs through vial adapter systems.

To find out how West's innovative systems can add value to your drug product or application, please contact your local West account representative.



Rx Only

The MixJect® Transfer Device is 510(k) cleared by the United States Food and Drug Administration and carries the CE mark (0344). Products are shown for INFORMATION purposes only and may not be approved for marketing in specific regions. Distribution and use are subject to applicable regulatory approvals and requirements for medical devices. The MixJect® transfer device is configurable and may not be suitable for use with all drugs. Refer to drug manufacturer's labeling and use instructions for device configuration compatibility. The device may contain a pre-attached needle containing cobalt. Injection site pain, irritation, discomfort, bruising, and/or erythema may occur. Failure to follow product instructions for use may result in compromised sterility; contamination; leakage (including possible exposure to medication); and/or inadequate medication reconstitution, transfer, and/or dosing. Product misuse could potentially lead to needlestick injury, user and/or patient exposure to pathogens or infection, and/or suboptimal or delayed therapy. Please contact your West Pharmaceutical Services, inc. (West) representative for product availability. Important product and safety information and warnings at: <https://www.westpharma.com/-/media/westpharma/files/products/mixject-indication-safety-and-warnings.pdf>

FOR MORE INFORMATION VISIT [WESTPHARMA.COM](https://www.westpharma.com)