



MIX2VIAL[®]

Needleless Transfer Device



The Mix2Vial[®] needleless transfer device enables simple, vial-to-vial transfer for the reconstitution of lyophilized drug products. The reconstituted drug is available for immediate aspiration into the syringe used for injection.

MIX2VIAL[®] NEEDLELESS TRANSFER DEVICE ATTRIBUTES

- Needle-free for safety
- Touch-free packaging maintains sterility
- Easy to use



PERFORMANCE OPTIONS

- In-line liquid filter
- Diluent flow restrictor to prevent foaming

The Mix2Vial[®] needleless transfer device is ideal for lyophilized drugs provided with fixed-volume diluent vials.

Components are available to work with 13mm, 17mm and 20mm stoppers.

INNOVATIONS IN SAFETY AND ADMINISTRATION SYSTEMS

West's administration systems enhance the safe delivery of drugs through vial adapter systems.

To find out how West's innovative systems can add value to your drug product or application, please contact your local West account representative.

The Mix2Vial[®] needleless transfer device is 510(k) cleared by the United States Food and Drug Administration and carries the CE mark (0344). Products are shown for INFORMATION purposes only and may not be approved for marketing in specific regions. Distribution and use are subject to applicable regulatory approvals and requirements for medical devices. Please contact your West Pharmaceutical Services, inc. (West) representative for product availability. Important product information, safety and warnings are included at: <https://www.westpharma.com/products/vial-adapter-systems/mix2vial-and-needle-free-transfer-device>

FOR MORE INFORMATION VISIT **WESTPHARMA.COM**

NORTH AMERICA 1-800-345-9800 OPTION 8 | **SOUTH AMERICA** +55 11 4055 6061

EUROPE +49 (0) 2403-7960 | **ASIA PACIFIC** +65 6862 3400

WEST PHARMACEUTICAL SERVICES, INC. 530 HERMAN O. WEST DR. | EXTON, PA 19341